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System Message	Type	Possible Cause	Corrective Action
Timed out while waiting for large enough buffer to resume consuming!	Warning	Instrument disconnected / power off during a run.	Reconnect / restart instrument
An error occurred during the run! The run will now be aborted	Warning	Instrument disconnected / power off during a run	Reconnect / restart instrument
Experiment New Experiment contains no acquisitions	Warning	The run does not contain acquisitions	Select an analysis mode and an acquisition type.
Please activate an instrument before setting up a new run	Warning	There is no default instrument defined.	Click on "Open Tools", goto Instruments and set an instrument as default.
Dispatcher Exception: Attempted to connect instrument xxxxx to a different instrument (yyyyy) on port HTC on 192.168.95.41. Please choose the correct instrument in the software, connect to a different instrument/ port, or create a new instrument if this is a new instrument (hasn't been connected before)!	Warning	The instrument connected to the IP address has changed.	Click on "Open Tools", goto Instruments and create a new instrument, and set it as default.
Failed to open log file C:\Program Files\Roche\LightCycler480\Bin\Logs\Instr_XXXXXX000.log. A connected instrument is already using the log file!	Warning	The instrument is already connected	-
Cannot find a user with the login name "xy"	Warning	User does not exist in database.	Click on "Open Tools" and check user settings under "Users and Groups". Define user as appropriate.
The supplied password does not match user xy!	Warning	Wrong password entered.	Check password entered!

System Message	Type	Possible Cause	Corrective Action
The application was unable to abort the run in a timely manner, please save the experiment to recover any data acquired this far. The application may not behave as expected until restarted!	Error	Hardware or software error in the instrument.	Save the experiment and restart LightCycler® 480 Software.
An error occurred during the run! The run will now be aborted	Error	Hardware or software error in the instrument.	If instrument's status LED is RED, restart instrument.

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Error No	System Message	Type	Possible Cause	Corrective Action
16	Restart instrument, inform Roche service	Error	Communication Error with Block Cycler Controller.	Restart instrument.
26	Restart instrument, inform Roche service	Error	Communication Error with Detection Controller.	Restart instrument.
36	Restart instrument, inform Roche service	Error	Communication Error with Raw Data controller 3.	Restart instrument.
46	Restart instrument, inform Roche service	Error	Communication Error with Detection Controller 3.	Restart instrument.
11001	Command not allowed	Error	The command send to the instrument cannot be executed at this time.	a) check MWP, try command again b) after error: restart LightCycler® Software and recover data from instrument c) restart instrument
10006	Wrong instrument type connected	Error	The software is not connected to a Light-Cycler® 480 Instrument or the instrument EEPROM is defect.	Check if you are connected to the correct instrument. Contact Roche technical service if it is the case.
100100010	Restart instrument, inform Roche service	Error	Block Cycler Controller Board removed or not properly connected.	Check Block Cycler Controller Board.
100100020	Restart instrument, inform Roche service	Error	Detection Control Motors Controller Board removed or not properly connected.	Check Detection Control Motors Controller Board.
100100040	Restart instrument, inform Roche service	Error	Detection Control Controller Board removed or not properly connected.	Check Detection Control Controller Board.

Error No	System Message	Type	Possible Cause	Corrective Action
100100050	Restart instrument, inform Roche service	Error	Data Reduction Controller Board removed or not properly connected.	Check Data Reduction Controller Board.
100100080	Restart instrument, inform Roche service	Error	Controller Firmware is older as expected.	Restart instrument.
100100100	Restart instrument, inform Roche service	Error	Indicator Board defect or unplugged.	Check Indicator Board.
100600120	Command not allowed	Warning	Command cannot be executed due to a previous error.	If instrument status LED is RED, restart instrument. If instrument status LED is GREEN, check multiwell plate and repeat command.
100600130	Command not allowed	Warning	Command cannot be executed due to a previous error.	If instrument status LED is RED, restart instrument. If instrument status LED is GREEN, check multiwell plate and repeat command.
100800100		Warning Error	Block cyler and cover are not compatible. Block cyler EEPROM defect.	Check thermal block cyler and block cyler cover. Check EEPROM cable.
101100100		Warning	Wrong MWP orientation. MWP was loaded during initialization.	Check multiwell plate.
101100110		Warning	Multiwell plate and block cyler are not compatible.	Check multiwell plate and thermal block cyler type.
101400130		Error	Door is open. Micro switch defect.	Close door. Check Micro switch.

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Error No.	System Message	Type	Possible Cause	Corrective Action
210259098	Restart Instrument, inform Roche service	Synch-Error	Block temperature out of range.	Check block cycler. Check settings for ramp rate (°C/s).
210259120	Restart Instrument, inform Roche service	Synch-Error	Electronic rack temperature over 90°C.	
210260004	Restart Instrument, inform Roche service	Synch-Error	Cover Heater on failed.	Check block cycler cover.
230257122	Restart Instrument, inform Roche service	Synch-Error	LB Ex Filter Wheel defect	Check Filter Wheel.
230258124	Restart Instrument, inform Roche service	Synch-Error	Motor Ex Filter Wheel defect	Check Filter Wheel.
230513130	Restart Instrument, inform Roche service	Synch-Error	Xenon Lamp defect.	Check Xenon lamp.
240117060	Restart Instrument, inform Roche service	Synch-Error	Camera temperature too high.	Check camera. Check environmental temperature.
310000096	Restart Instrument, inform Roche service	Asynch-Error	No block cycler inserted. Sensor data are not available.	Check Block, Block Connection.
310000099	Restart Instrument, inform Roche service	Asynch-Error	No block cycler cover or not recognized.	Check block cycler cover heater.
310000107	Restart Instrument, inform Roche service	Asynch-Error		Check block cycler.
310000108		Asynch-Error	Block cycler temperature sensor defect.	Check block cycler.
310000112	Restart Instrument, inform Roche service	Asynch-Error	No block cycler inserted. Sensor data are not available.	Check block cycler and block connection.
310000115	Restart Instrument, inform Roche service	Asynch-Error	No block inserted Sensor data are not available.	Check block cycler and block connection.
310000118	Restart Instrument, inform Roche service	Asynch-Error	Electronic Rack temperature sensor unplugged or defect	Check Electronic Rack temperature sensor.
310000145	Restart Instrument, inform Roche service	Asynch-Error	Fuse F5 defect.	Check fuse F5.
310000146	Restart Instrument, inform Roche service	Asynch-Error	Fuse F4 defect.	Check fuse F4.
310000147	Restart Instrument, inform Roche service	Asynch-Error	Fuse F3 defect.	Check fuse F3.
310000148	Restart Instrument, inform Roche service	Asynch-Error	Fuse F2 defect.	Check fuse F2.

Error No.	System Message	Type	Possible Cause	Corrective Action
340000050	Restart Instrument, inform Roche service	Asynch-Error	Camera FireWire cable unplugged. Camera Datacable unplugged. Ex Motor or LB on Optics unit defect. Electronic Rack Fan defect.	Check Camera. Check Optics module. Check Electronic Rack Fan.
340000060	Restart Instrument, inform Roche service	Asynch-Error	Camera temperature too high.	Check camera. Check environmental temperature.
340000067	Restart Instrument, inform Roche service	Asynch-Error	Cover Markers not found.	Check block cyclor cover and camera.
340000080	Restart Instrument, inform Roche service	Asynch-Error	Optics cable unplugged. LB Em Filter Wheel defect. Block cyclor cover markers not found. The Reference Channel EEPROM is not programmed.	Check Xenon lamp. Check Optics module. Check Filter Parameters.
500100850	Restart Instrument, inform Roche service	Error	Appears in combination with 230513130 (Xenon Lamp defect).	Check Xenon lamp.
500500420	Restart Instrument, inform Roche service	Error	No open channel Found TL-Slave 1.	Check instrument connection.

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